



Joint Workshop of Junior Research - JWJR



Environmental, social and economic sustainability

EVALUATION OF AN ENVIRONMENTAL FRIENDLY METHOD FOR THE POLYMERIC E-WASTE DIGESTION FOR FURTHER Br AND METALS DETERMINATION

Gabrielle D. Iop (PG), Sindy R. Krzyzaniak (PG), Alice P. Holkem (PG), Gustavo R. Bitencourt (PG),
Jussiane S. Silva (PQ), Érico M. M. Flores (PQ), and Paola A. Mello (PQ)

Universidade Federal de Santa Maria, Departamento de Química, Santa Maria, RS - Brazil


May 10th, 2019

ELECTRICAL AND ELECTRONIC EQUIPMENT



Essential to our modern way of life

➤ **Mobile phone, computer, TV, among others.**



- **Improved performance and reduced cost in each new generation of product**
- **Encourage unsustainable behaviour**
- **Mobile phones: often treated as fashion items replaced long before their design lifetimes have expired**

Electronic Waste (E-Waste)



ELECTRONIC WASTE

- Discharge
- Recycling



Waste management is difficult:

- Hazardous components
- Formation of hazardous substances

**Legislation to overcome the increase
in generated e-waste**





- **Brazil: there is no specific policy for this type of material**
 - **CONAMA: disposal of batteries**
 - **National Policy on Solid Waste Management: producer responsibility**

ELECTRONIC WASTE

Composition

- ❖ Printed circuits boards
- ❖ Polymers or blends
- ❖ Additives (chemical, thermic and mechanical resistance)
- ❖ Brominated compounds (polybrominated biphenyls - PBBs- and polybrominated diphenylethers -PBDEs-)



- **Additives: contain toxic elements, as Sb, Cd, Hg**
- ↪ **Toxic to the environment and human lives**

ELECTRONIC WASTE



Legislation: Europe

❖ Directive 2011/65/EU (Restriction of the Use of Certain Hazardous Substances, RoHS)

- Cd (0.01%)
- Cr(VI), Hg, and Pb (0.1%)
- Polybrominated compounds (0.1%)

→ PBB and PBDE


➤ Sb: not specified in the Directive
→ Used as flame retardant

➤ Polybrominated compounds: about 50% (m/m) of Br
→ Bromine determination
(indirect determination of compounds)


BROMINE AND METALS DETERMINATION IN POLYMERS



SAMPLE IN SOLUTION



SAMPLE PREPARATION STEP



- ❖ **HIGH STABILITY**
- ❖ **RESISTANCE TO ACID ATTACK**
- ❖ **DIGESTS MUST BE SUITABLE FOR THE DETERMINATION TECHNIQUES**

↪ **Use of concentrated acids**



SAMPLE IN
SOLU

SAMPLE
ATION STEP

Analytical Challenge



- ❖ HIGH STABILITY
- ❖ RESISTANCE TO ACID ATTACK
- ❖ DIGESTS MUST BE SUITABLE FOR THE DETERMINATION TECHNIQUES

↪ Use of concentrated acids

Sample preparation step

WET DIGESTION

EXTRACTION

MIC

- ✓ Fast heating of the solution
- ✓ High temperature
- ✗ Concentrated acids
- ✗ High carbon content
- ✗ High residual acidity

- ✓ High sample mass
- ✓ High temperature
- ✗ Reproducibility
- ✗ Extraction efficiency
- ✗ High carbon content


- ✓ Low carbon content
- ✓ Diluted solutions
- ✗ Cost

Problem: there is no single protocol for bromine and metals

MIC: Microwave-Induced Combustion

SAMPLE PRETREATMENT AND SEPARATION

Greener sample preparation procedures




MIC: Microwave-Induced Combustion

ICP: Inductively Coupled Plasma

UV: Ultraviolet

SAMPLE PRETREATMENT AND SEPARATION

Greener sample preparation procedures




MIC: Microwave-Induced Combustion

ICP: Inductively Coupled Plasma

UV: Ultraviolet

OBJECTIVE

- To develop a single sample preparation method for polymeric e-waste digestion by microwave-assisted ultraviolet digestion (MAWD-UV)
- To use of diluted acids for further bromine and metals determination by inductively coupled plasma mass spectrometry (ICP-MS)



ELECTRONIC WASTE (E-Waste)

Sample preparation methods

MAWD-UV

Sample mass: 200 mg

Digestion mixture (15 mL):

- A. HNO_3 (1 to 7 mol L^{-1})
- B. HNO_3 + HCl (6 + 1, 3 + 0.5, 1.2 + 0.2 or 0.6 + 0.1 mol L^{-1})
- C. 1 mol L^{-1} HNO_3 + H_2O_2 (0, 1.5, 3.5 or 4.7 mol L^{-1})
- D. 0.5 mol L^{-1} HNO_3 + H_2O_2 (0, 1.5, 3.5, 4.7 or 5.3 mol L^{-1})
- E. HNO_3 + HCl (0 + 0.5, 0.5 + 0.5, 0.5 + 1 or 0.5 + 0 mol L^{-1}) + 5.3 mol L^{-1} H_2O_2

Heating program:

- i) 700 W for 40 min (ramp of 25 min)
Exhaustion at 125 m^3 air h^{-1}
 - ii) 0 W for 20 min for cooling
Exhaustion at 125 m^3 air h^{-1}
- Max. pressure: 80 bar
Max. temperature: 250 °C

Dilution (25 mL)

Determination of Br, Cd, Cr, Pb, and Sb by ICP-MS and Hg by FI-CVG-ICP-MS

MAWD

Sample mass: 200 mg

Digestion mixture (15 mL):

- A. HNO_3 (1 to 7 mol L^{-1})
- B. HNO_3 + HCl (6 + 1, 3 + 0.5, 1.2 + 0.2 or 0.6 + 0.1 mol L^{-1})

Heating program:

- i) 700 W for 40 min (ramp of 25 min)
Exhaustion at 60 m^3 air h^{-1}
 - ii) 0 W for 20 min for cooling
Exhaustion at 125 m^3 air h^{-1}
- Max. pressure: 80 bar
Max. temperature: 280 °C

Dilution (25 mL)

Reference values

NAA for Br, Cr, and Sb

MAWD for Cd, Hg, and Pb

Sample mass: 200 mg

Digestion mixture (15 mL):

- 5 mL 14 mol L^{-1} HNO_3 + 1 mL 12 mol L^{-1} HCl

Heating program:

- i) 1400 W for 35 min (ramp of 10 min)
Exhaustion at 60 m^3 air h^{-1}
 - ii) 0 W for 20 min for cooling
Exhaustion at 125 m^3 air h^{-1}
- Max. pressure: 80 bar
Max. temperature: 280 °C

Determination of Cd and Pb by ICP-MS and Hg by FI-CVG-ICP-MS

MAWD-UV: Microwave-Assisted Ultraviolet Digestion

MAWD: Microwave-assisted Wet Digestion


NAA: Neutron Activation Analysis

ICP-MS: Inductively Coupled Plasma Mass Spectrometry

FI-CVG-ICP-MS: Flow-Injection Cold Vapor Generation Inductively Coupled Plasma Mass Spectrometry


EXPERIMENTAL

Grinding polymeric parts with cryogenic mill




Cryogenic mil:
2 min of pre-cooling followed by 3 min for grinding
(last procedure repeated three times)

Microwave-assisted ultraviolet digestion (MAWD-UV)



UV lamp and PTFE devices
1- lamp base ring
2- lamp spacer



Lamp base ring, 5
mL of solution and
200 mg of sample



UV lamp with the
lamp spacer

10 mL of
solution (total
volume: 15 mL)

MW irradiation

RESULTS


➤ Influence of UV radiation for e-waste digestion



■ MAWD
■ MAWD-UV

● MAWD
■ MAWD-UV



Improvement in digestion efficiency when more concentrated HNO₃ solutions and more concentrated HNO₃ + HCl mixtures were used



Digestion efficiency was more pronounced when the UV lamp was used

RESULTS

- Influence of the digestion solution in MAWD-UV
 - 1 to 7 mol L⁻¹ HNO₃ and mixtures of HNO₃ + HCl



All conditions resulted in lower values for bromine in comparison with the reference value



UV lamp: improved the destruction of the organic matrix

RESULTS

➤ Influence of the digestion solution in MAWD-UV

- Acid mixture with $5.3 \text{ mol L}^{-1} \text{ H}_2\text{O}_2$

Analyte	$\text{HNO}_3 + \text{HCl}, \text{ mol L}^{-1}$				Reference value ^{a,b}
	0 + 0.5	0.5 + 0	0.5 + 0.5	0.5 + 1	
Br, mg g^{-1}	1.72 ± 0.19	8.32 ± 0.68	8.82 ± 0.67	5.76 ± 0.04	8.8 ± 0.2^a
Cd, $\mu\text{g g}^{-1}$	8.16 ± 0.31	14.4 ± 1.2	16.0 ± 1.0	10.9 ± 0.1	15.8 ± 0.3^b
Cr, $\mu\text{g g}^{-1}$	3.32 ± 0.37	3.02 ± 0.55	6.30 ± 0.70	5.45 ± 0.33	9 ± 1^a
Hg ^c , ng g^{-1}	62.0 ± 3.0	15.0 ± 1.0	66.2 ± 2.6	50.7 ± 4.3	65.9 ± 4.5^b
Pb, $\mu\text{g g}^{-1}$	11.5 ± 0.8	25.1 ± 1.8	24.9 ± 1.4	21.4 ± 1.1	27.3 ± 2.6^b
Sb, mg g^{-1}	2.17 ± 0.07	2.62 ± 0.26	2.99 ± 0.24	3.10 ± 0.32	3.27 ± 0.09^a


^a Results obtained by NAA.

^b Results obtained after MAWD.

^c Determination by FI-CVG-ICP-MS.

Considered as a compromise condition to provide quantitative recoveries

- Carbon content lower than 300 mg L^{-1}
- Residual acidity about 70%



Cr: no acid mixture was suitable for its quantitative recovery

Form of Cr present in the sample is important:

Use of extreme conditions for digestion (temperature and pressure)

RESULTS

Accuracy: CRM of low-density polyethylene (ERM 680k and ERM 681k)

Analyte	ERM EC680k		ERM EC681k	
	Found value	Certified value	Found value	Certified value
Br, $\mu\text{g g}^{-1}$	97 ± 8	96 ± 4	725 ± 71	770 ± 40
Cd, $\mu\text{g g}^{-1}$	18.3 ± 0.5	19.6 ± 1.4	132 ± 9	137 ± 4
Cr, $\mu\text{g g}^{-1}$	4.14 ± 0.27	20.2 ± 1.1	24.5 ± 6.1	100 ± 5
Hg*, $\mu\text{g g}^{-1}$	4.63 ± 0.22	4.64 ± 0.20	23.3 ± 0.6	23.7 ± 0.8
Pb, $\mu\text{g g}^{-1}$	12.5 ± 0.9	13.6 ± 0.5	91.2 ± 6.1	98 ± 6
Sb, $\mu\text{g g}^{-1}$	9.59 ± 0.19	10.1 ± 1.6	97.4 ± 8.0	99 ± 6


* Determination by FI-CVG-ICP-MS

Digestion of 200 mg of sample, n = 3.

Results for all elements by MAWD-UV were obtained using 15 mL of $0.5 \text{ mol L}^{-1} \text{ HNO}_3 + 0.5 \text{ mol L}^{-1} \text{ HCl} + 5.3 \text{ mol L}^{-1} \text{ H}_2\text{O}_2$ for digestion.

**No statistical difference in comparison with the certified values
(*t*-test, 95% confidence level)**

Spike recoveries



From 90% (for Hg) to 113% (for Sb)

CONCLUSION

- ❖ The developed method can be considered **environmental friendly**: use of diluted acid mixture, reduced time, less reagents consumption, and energy saving
- ❖ Possibility of using a **single protocol** for further Br and metals determination

ACKNOWLEDGEMENTS

Organizing Committee



Asociación de Universidades
GRUPO MONTEVIDEO



Funding Agencies



SAMPLE PREPARATION STEP

Digestion efficiency

- ✓ **Residual acidity**

% , in relation to the initial amount of acid

- ✓ **Residual carbon content (RCC)**

% , in relation to the amount of C originally present in the sample


- ✓ **Carbon content in digests/Dissolved organic carbon (DOC)**

mg L⁻¹, dissolved carbon in final digests

Suitable digests: ↓ RA, RCC and DOC

Procedures for e-waste digestion

Sample preparation step



Problem: there is no single protocol for bromine and metals

MIC: Microwave-Induced Combustion